

# Actimove<sup>®</sup>

GenuFlex



- Non slip, wrap-around design
- Comfortable, breathable fabric
- Advanced polycentric hinges

## Functional hinged knee brace

Feel good. Move better.

# Actimove® GenuFlex

## A balanced concept.

Actimove® GenuFlex is constructed from advanced breathable, neoprene-free fabric which provides effective compression and pain relief. The slim-line, wrap-around design makes it easy to apply and wear under normal clothes whilst preventing brace migration. Knee movements are guided anatomically by polycentric hinges, providing reassuringly safe knee support.

### Indications

- Minor ligament instabilities of the knee
- Chronic meniscal damage
- Post-op use
- Painful osteoarthritis of the knee
- Osteoarthritis and rheumatoid arthritis

### Stabilisation

- Neoprene-free AirFlow™ material provides good compression
- Polycentric hinges guide knee movements anatomically allowing full flexion and preventing hyperextension
- Advanced counter-strap system reduces brace migration and rotation

### Comfort


- Breathable, moisture transmitting AirFlow™ material
- Snug, non-slip fit for enduring comfort and strong support
- Special closure system allows quick and easy application
- Open popliteal area prevents wrinkling and skin maceration
- Molded hinges prevent pressure points
- Can be worn under everyday clothes due to its slim design

### Design

- Wrap-around design for easy, convenient application
- Patella opening helps correct positioning
- Attractive modern design enhances patient compliance



### Ordering information



BSN Code	NHS Code	Size	Thigh Circumference*		Sales Unit
73492-32	GRD059	XS	30–34 cm	12–13 1/2"	1/1
73492-33	GRD060	S	35–39 cm	14–15 1/2"	1/1
73492-34	GRD061	M	40–44 cm	16–17 1/2"	1/1
73492-35	GRD062	L	45–49 cm	18–19 1/2"	1/1
73492-36	GRD063	XL	50–54 cm	20–21 1/2"	1/1
73492-37	GRD064	XXL	55–59 cm	22–23 1/2"	1/1
73492-38	GRD065	XXXL	60–64 cm	24–25 1/2"	1/1

\*Circumference measured 10 cm (4 inch) above the kneecap with a flexion of 15°



**LATEX  
FREE  
FORMULA**



BSN medical Ltd, PO BOX 258, Willerby, Hull HU10 6WT  
Tel. 01482 670100, Fax: 01482 670111, Email: orders.uk@bsnmedical.com www.bsnmedical.com